

Product datasheet for **MP210527**

Kat2b Mouse qPCR Primer Pair (NM_020005)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18519
Forward Sequence:	CCTCTTCACCTGCGTCCACAAA
Reverse Sequence:	TCTCCAAGGAGCCTTCAACCAC
ACCN:	BC082581 , NM_020005 , NM_020005.1 , NM_020005.2 , NM_020005.3 , NM_020005.4 , BC019722 , BC138195 , BC145896
UniProt ID:	Q9JHD1
Synonyms:	A930006P13Rik; AI461839; AW536563; p/CAF; Pcaf
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)