

Product datasheet for **MP210531**

Paf1 Mouse qPCR Primer Pair (NM_019458)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54624
Forward Sequence:	GGAGGAAGAGATGGAGGCTGAA
Reverse Sequence:	CACTTGCCCTCATCTCTGTCACC
ACCN:	<u>BC029843</u> , <u>BC083337</u> , <u>NM_019458</u> , <u>NM_019458.1</u> , <u>NM_019458.2</u> , <u>NM_019458.3</u> , <u>BC010317</u>
UniProt ID:	<u>Q8K2T8</u>
Synonyms:	5730511K23Rik; AA675045
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)