

Product datasheet for **MP210559**

Pak1ip1 Mouse qPCR Primer Pair (NM_026550)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68083
Forward Sequence:	GGATGTGAAGAGATGGAAGTGCC
Reverse Sequence:	TATCGGTGCCAACAGACAGAGC
ACCN:	<u>NM_026550</u> , <u>NM_026550.1</u> , <u>NM_026550.2</u> , <u>NM_026550.3</u> , <u>BC018522</u> , <u>BC024629</u> , <u>BC043461</u> , <u>BC043725</u> , <u>BC051498</u> , <u>BC060371</u> , <u>BC116729</u> , <u>BC116731</u> , <u>BY055321</u> , <u>BY788898</u>
UniProt ID:	<u>Q9DCE5</u>
Synonyms:	5830431I15Rik; 5930415H02Rik; AA419825; AI314040; AW556169; Gdpd1; PIP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)