

## Product datasheet for **MP210580**

### Pank4 Mouse qPCR Primer Pair (NM\_172990)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	269614
Forward Sequence:	GCTCGTGGATTCTACACCAAG
Reverse Sequence:	CCCTGACAAAGGTGAAGACTCC
ACCN:	<u><a href="#">BC050089</a></u> , <u><a href="#">NM_172990</a></u> , <u><a href="#">NM_172990.1</a></u> , <u><a href="#">NM_172990.2</a></u> , <u><a href="#">NM_172990.3</a></u> , <u><a href="#">NM_172990.4</a></u> , <u><a href="#">NM_172990.5</a></u> , <u><a href="#">BC050089.1</a></u> , <u><a href="#">BC012520</a></u> , <u><a href="#">BC019659</a></u> , <u><a href="#">BC027074</a></u> , <u><a href="#">BK000016</a></u>
UniProt ID:	<u><a href="#">Q80YV4</a></u>
Synonyms:	D030031112Rik; R75150
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)