

Product datasheet for **MP210586**

Pard3b Mouse qPCR Primer Pair (NM_001081050)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72823
Forward Sequence:	CCAGTTGGTGATTCTGAGGAGAC
Reverse Sequence:	ATTTCGCCTGGAGGTGCTCAGA
ACCN:	<u>NM_001081050</u> , <u>NM_001081050.1</u> , <u>NM_001081050.2</u> , <u>BC008193</u> , <u>BC114605</u> , <u>BC167176</u> , <u>BF461790</u>
UniProt ID:	<u>Q9CSB4</u>
Synonyms:	1810008K04Rik; 2010002N16Rik; 2810455B10Rik; Als2cr19; PAR3B; PAR3beta; PAR3L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)