

Product datasheet for **MP210588**

Pard6a Mouse qPCR Primer Pair (NM_019695)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56513
Forward Sequence:	TACGGATGCTCACGGTGACTTG
Reverse Sequence:	GAGAGTTGGAGGCAAAAGCCAG
ACCN:	BC049593 , BC061114 , NM_019695 , NM_019695.1 , NM_019695.2 , NM_019695.3 , BF466814 , BG806503 , BY099545 , BY280124
UniProt ID:	Q9Z101
Synonyms:	0710008C04Rik; 2610010A15Rik; Par-6; PAR-6A; Par6; PAR6alpha; Par6c; TAX40; Tip-40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)