

Product datasheet for **MP210612**

Palb2 Mouse qPCR Primer Pair (NM_001081238)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	233826
Forward Sequence:	CATACCTGGCACTTCACAGAGG
Reverse Sequence:	GTCATCACCAGAAGAGCACAGC
ACCN:	BC066140 , NM_001081238 , NM_001081238.1 , NM_001081238.2 , BB061170 , BB757034 , BC044748 , BC055302 , BC066139 , BQ560264
UniProt ID:	Q3U0P1
Synonyms:	BC066140
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)