

Product datasheet for **MP210624**

Pnpla7 Mouse qPCR Primer Pair (NM_146251)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	241274
Forward Sequence:	ACCATCAGAGCCAACCGAGACT
Reverse Sequence:	CTGCCGTACAAAGGATGACATCC
ACCN:	<u>NM_146251</u> , <u>NM_146251.1</u> , <u>NM_146251.2</u> , <u>NM_146251.3</u> , <u>NM_146251.4</u> , <u>BC141365</u> , <u>BB618610</u> , <u>BB852360</u> , <u>BC021312</u> , <u>BC025621</u> , <u>BC027342</u> , <u>BC064003</u> , <u>BC080793</u> , <u>BC141366</u> , <u>BE623320</u>
UniProt ID:	<u>A2AJ88</u>
Synonyms:	BC027342; E430013P11Rik; Nre
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)