

Product datasheet for **MP210639**

Cdk17 Mouse qPCR Primer Pair (NM_146239)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237459
Forward Sequence:	GTCACCTACAGGTGTTTGCCTC
Reverse Sequence:	TGCAGCTTCTCCAGGTAGCCAT
ACCN:	BC031778 , BC049904 , BC064815 , NM_146239 , NM_146239.1 , NM_146239.2 , BB085519 , BQ442862 , BU841444 , BY130550
UniProt ID:	Q8K0D0
Synonyms:	6430598J10Rik; Pctk2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)