

Product datasheet for **MP210649**

Pdlim5 Mouse qPCR Primer Pair (NM_019808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56376
Forward Sequence:	TGTGTTCCGAGTCTGCTCAGGA
Reverse Sequence:	ACTCCGTGTTGCGCTCCACAAT
ACCN:	BC037476 , NM_019808 , NM_019808.1 , NM_019808.2 , NM_019808.3 , BC020145 , BG064526 , BG077941 , BY227455 , BY741234
UniProt ID:	Q8CI51
Synonyms:	1110001A05Rik; AI987914; C87059; Enh; Enh1; Enh2; Enh3; LIM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)