

Product datasheet for **MP210661**

Pcnx Mouse qPCR Primer Pair (NM_018814)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54604
Forward Sequence:	GAGCACAGTGTTACGGCATTG
Reverse Sequence:	TTCCTCACCTGAGGCAGCACAT
ACCN:	<u>NM_018814</u> , <u>NM_018814.1</u> , <u>NM_018814.2</u> , <u>NM_018814.3</u> , <u>BC028455</u> , <u>BC057333</u> , <u>BC066849</u> , <u>BC094637</u> , <u>BC098370</u> , <u>BC156548</u> , <u>BC172701</u> , <u>BP753289</u>
UniProt ID:	<u>Q9QYC1</u>
Synonyms:	2900024E21Rik; 3526401J03Rik; AF096286; AI327143; AI413187; Pcnx1; Pcnxl1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)