

Product datasheet for **MP210666**

Peli3 Mouse qPCR Primer Pair (NM_172835)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	240518
Forward Sequence:	CTGCTGGCTTTGATGCCTCTAG
Reverse Sequence:	GCAGAGTCCTCGGAGAAGCCA
ACCN:	BC028931 , NM_172835 , NM_172835.1 , NM_172835.2 , NM_172835.3 , BC022735 , BC057046 , BY128359 , BY452323
UniProt ID:	Q8BXR6
Synonyms:	6030441F14Rik; A930011L17; BC028931
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)