

Product datasheet for **MP210668**

Pef1 Mouse qPCR Primer Pair (NM_026441)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67898
Forward Sequence:	GGAAATTCCTCCAGCAGTGGAG
Reverse Sequence:	AACTGAGGGCTCAGGTTGTAGC
ACCN:	<u>BC019191, NM_026441, NM_026441.1, NM_026441.2, NM_026441.3, NM_026441.4</u>
UniProt ID:	<u>Q8BFY6</u>
Synonyms:	2600002E23Rik; Peflin
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)