

Product datasheet for **MP210671**

Ptcd3 Mouse qPCR Primer Pair (NM_027275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69956
Forward Sequence:	CATCCAAGATGGAGTCTGGTCC
Reverse Sequence:	CTCGGATCATCGTGCAGTAGGA
ACCN:	BC115434 , NM_027275 , NM_027275.1 , NM_027275.2 , NM_027275.3 , BB844000 , BC026759 , BC055844 , BC060687 , BC068263 , BC084595 , BF020753
UniProt ID:	Q14C51
Synonyms:	2610034F17Rik; 2810422B04Rik; AA589622; AU045708; MRP-S39
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)