

Product datasheet for **MP210674**

Pmpca Mouse qPCR Primer Pair (NM_173180)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66865
Forward Sequence:	CTGTTGGCTGATGTGGTTCTGC
Reverse Sequence:	GAATCATCTCGGTGAGCAGTGG
ACCN:	<u>BC010810</u> , <u>NM_173180</u> , <u>NM_173180.1</u> , <u>NM_173180.2</u> , <u>NM_173180.3</u>
UniProt ID:	<u>Q9DC61</u>
Synonyms:	I200002L24Rik; 4933435E07Rik; Alpha-MPP; P-55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

