

Product datasheet for **MP210695**

Ppil3 Mouse qPCR Primer Pair (NM_027351)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70225
Forward Sequence:	GAAGAGGAGGTAGCAGCATCTG
Reverse Sequence:	TGTGGCTGCTTGCCATAGGTGA
ACCN:	<u>BC060993</u> , <u>BC061645</u> , <u>NM_027351</u> , <u>NM_027351.1</u> , <u>NM_027351.2</u> , <u>NM_027351.3</u> , <u>BB687332</u> , <u>BB844274</u> , <u>BY288296</u>
UniProt ID:	<u>Q9D6L8</u>
Synonyms:	2310076N22Rik; 2510026K04Rik; Cyp10l
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)