

Product datasheet for **MP210704**

Ppig Mouse qPCR Primer Pair (NM_001081086)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	228005
Forward Sequence:	CCTCCTAAAGCCGATGACAAGG
Reverse Sequence:	ATCGTCGCTGGTATGAAGCAGG
ACCN:	NM_001081086 , NM_001081086.1 , BC150694 , BC016126 , BC027420 , BC038818 , BC048905 , BC058789 , BC094926
UniProt ID:	A2AR02
Synonyms:	AU019516; AU022200; B230312B02Rik; CYP; SRCyp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)