

Product datasheet for **MP210714**

Pcnt Mouse qPCR Primer Pair (NM_008787)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18541
Forward Sequence:	GAGGAGAAGTCGGTCTTGTGGA
Reverse Sequence:	GCGGTCCTTTTCAGACTGCTTC
ACCN:	<u>NM_008787</u> , <u>NM_008787.1</u> , <u>NM_008787.2</u> , <u>NM_008787.3</u> , <u>BB840023</u> , <u>BC034865</u> , <u>BC039661</u> , <u>BC055467</u> , <u>BC085244</u> , <u>BC100380</u> , <u>BC172748</u> , <u>BG084686</u> , <u>BG088995</u> , <u>BG143796</u> , <u>BG923243</u> , <u>BI734431</u> , <u>BM219062</u> , <u>BU518252</u> , <u>BY711114</u>
UniProt ID:	<u>P48725</u>
Synonyms:	AW476095; C86676; KEN; kendrin; m239Asp; m275Asp; Pcnt2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)