

Product datasheet for **MP210720**

Postn Mouse qPCR Primer Pair (NM_015784)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50706
Forward Sequence:	CAGCAAACCACTTTCACCGACC
Reverse Sequence:	AGAAGGCGTTGGTCCATGCTCA
ACCN:	BC031449 , NM_015784 , NM_015784.1 , NM_015784.2 , NM_015784.3 , BC007141 , BE850247 , BF020820
UniProt ID:	Q62009
Synonyms:	A630052E07Rik; AI747096; Os; OSF; OSF-2; Osf2; p; Pe; PLF; PN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)