

Product datasheet for **MP210751**

Pex16 Mouse qPCR Primer Pair (NM_145122)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18633
Forward Sequence:	CTGGCTACTGTCTGGTGTGGTA
Reverse Sequence:	AGAACGGTGAGCGTAGCAGGTA
Accession No:	<u>BC010822</u> , <u>NM_145122</u> , <u>NR_028114</u> , <u>BE690943</u>
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)