

Product datasheet for **MP210767**

Pcnp Mouse qPCR Primer Pair (NM_001024622)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76302
Forward Sequence:	GAAAGCGTCCGCCATCTCCATT
Reverse Sequence:	GGAGGCATTCTTCTGGCTCAC
ACCN:	BC061235 , BC132588 , BC132590 , NM_001024622 , NM_001024622.2 , BC023320 , BC106192 , BC145194 , BY287179
UniProt ID:	Q6P8I4
Synonyms:	1110018D06Rik; AI647035
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

