

Product datasheet for **MP210776**

Phf14 Mouse qPCR Primer Pair (NM_029404)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75725
Forward Sequence:	GTTCTGGAGTGGCTTTTGGAG
Reverse Sequence:	CCAGTTCTAGCAAAGCGAGGATC
ACCN:	<u>BC040236</u> , <u>NM_029404</u> , <u>NM_029404.1</u> , <u>NM_029404.2</u> , <u>NM_029404.3</u> , <u>BC040236.1</u> , <u>BC040235</u> , <u>BX524731</u> , <u>NM_029404.4</u>
UniProt ID:	<u>Q9D4H9</u>
Synonyms:	1110001C23Rik; 4932409F11Rik; 5730446A07Rik; AA623952; AV297001; mKIAA0783
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

