

Product datasheet for **MP210787**

Phf3 Mouse qPCR Primer Pair (NM_001081080)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	213109
Forward Sequence:	CAGAGAGAAGTGGAAAGACGACC
Reverse Sequence:	CTCAGACACATCTGGCTTCTCC
ACCN:	<u>NM_001081080</u> , <u>NM_001081080.1</u> , <u>BC137913</u> , <u>BC011226</u> , <u>BC019710</u> , <u>BC032205</u> , <u>BC038053</u> , <u>BC099524</u> , <u>BG069249</u>
Synonyms:	2310061N19Rik; mKIAA0244
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

