

Product datasheet for **MP210807**

Phactr4 Mouse qPCR Primer Pair (NM_175306)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100169
Forward Sequence:	AGCTTCTCCGAGGACACTGGTA
Reverse Sequence:	TAGGAGTGGCTGAGGTAGACGT
ACCN:	<u>BC075672</u> , <u>BC096033</u> , <u>NM_175306</u> , <u>NM_175306.1</u> , <u>NM_175306.2</u> , <u>NM_175306.3</u> , <u>NM_175306.4</u> , <u>BB849287</u> , <u>BC037059</u> , <u>BY294063</u>
UniProt ID:	<u>Q501J7</u>
Synonyms:	3110001B12Rik; AI527228; AW495572; C330013F19Rik; mKIAA4120; N28169
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

