

Product datasheet for **MP210808**

Pten Mouse qPCR Primer Pair (NM_008960)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19211
Forward Sequence:	TGAGTTCCTCAGCCATTGCCT
Reverse Sequence:	GAGGTTTCCTCTGGTCTGGTA
ACCN:	<u>BC021445</u> , <u>NM_008960</u> , <u>NM_008960.1</u> , <u>NM_008960.2</u> , <u>BC150968</u> , <u>BK008841</u> , <u>BU701098</u>
UniProt ID:	<u>O08586</u>
Synonyms:	2310035O07Rik; A130070J02Rik; A1463227; B430203M17Rik; MMAC; MMAC1; PTENbeta; TEP; TEP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)