

Product datasheet for **MP210816**

Plpp5 Mouse qPCR Primer Pair (NM_028000)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71910
Forward Sequence:	CGTGTGACTACAAGCACCCTG
Reverse Sequence:	GACGGAAGTTTGTGCTTGTCTG
ACCN:	<u>NM_028000</u> , <u>NM_028000.1</u> , <u>BB857717</u> , <u>BC023111</u> , <u>BC039626</u> , <u>BC062173</u> , <u>BC089556</u> , <u>BC099489</u> , <u>BE646843</u>
UniProt ID:	<u>Q3UMZ3</u>
Synonyms:	1810019D05Rik; 2310022A04Rik; AW045471; Ppapdcl; Ppapdclb
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

