

Product datasheet for **MP210818**

Pipp6 Mouse qPCR Primer Pair (NM_028922)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74411
Forward Sequence:	ATCGACCTGTGGCTGTCCAAGA
Reverse Sequence:	TCGCTCCTAAGCAGGCAGTAGA
ACCN:	BC052412 , NM_028922 , NM_028922.1 , NM_028922.2 , NM_028922.3 , BC051510
UniProt ID:	Q9D4F2
Synonyms:	4932443D16Rik; Ppapdc2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

