

Product datasheet for **MP210820**

Plpp3 Mouse qPCR Primer Pair (NM_080555)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67916
Forward Sequence:	CGGGATTGTCACGGGTATCTGA
Reverse Sequence:	GCAGTGAGAGGCTCGTCTTAGT
ACCN:	<u>BC005558</u> , <u>NM_080555</u> , <u>NM_080555.1</u> , <u>NM_080555.2</u> , <u>BY129827</u>
UniProt ID:	<u>Q99JY8</u>
Synonyms:	1110003O22Rik; 2610002D05Rik; AV025606; D4Bwg0538e; D4Bwg1535e; Lpp3; Ppab2b; Ppap2b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)