

## Product datasheet for **MP210821**

### Plpp2 Mouse qPCR Primer Pair (NM\_015817)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50784
Forward Sequence:	CATACCGTCCAGACACGATCAC
Reverse Sequence:	GGATCGCGAATAAAGACGGTCTG
ACCN:	<a href="#">BC010332</a> , <a href="#">NM_015817</a> , <a href="#">NM_015817.1</a> , <a href="#">NM_015817.2</a> , <a href="#">NM_015817.3</a> , <a href="#">BB737982</a> , <a href="#">BB868752</a> , <a href="#">BC016635</a> , <a href="#">BC085104</a> , <a href="#">BQ571983</a> , <a href="#">BY344892</a>
UniProt ID:	<a href="#">Q9DAX2</a>
Synonyms:	L; Lpp2; PAP2-G; PAP2-gamma; PAP2c; Pp; Ppap2c
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)