

Product datasheet for **MP210841**

Picalm Mouse qPCR Primer Pair (NM_146194)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 233489 |
| Forward Sequence: | GATGTTACCTGCCATTGCTTC |
| Reverse Sequence: | TGGCTGTGCAACTGGAGAAGGA |
| ACCN: | <u>BC011470</u> , <u>BC021491</u> , <u>BC023843</u> , <u>BC025566</u> , <u>BC057683</u> , <u>NM_146194</u> , <u>NM_146194.1</u> , <u>NM_146194.2</u> , <u>NM_146194.3</u> , <u>NM_146194.4</u> , <u>BC011470.1</u> , <u>BC025566.1</u> , <u>BC027116</u> , <u>BK001028</u> , <u>BY160265</u> , <u>NM_146194.5</u> |
| UniProt ID: | <u>Q7M6Y3</u> |
| Synonyms: | CALM; CLTH; fit-1; fit1; mKIAA4114; PAP180 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)