

Product datasheet for **MP210856**

Pigt Mouse qPCR Primer Pair (NM_133779)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78928
Forward Sequence:	ACGCTCACCATTACCTCCAAGG
Reverse Sequence:	ACTGGATGGAGACCTTGGTGAC
ACCN:	<u>BC034175</u> , <u>NM_133779</u> , <u>NM_133779.1</u> , <u>NM_133779.2</u> , <u>BC016520</u> , <u>BE994563</u> , <u>BU058874</u>
UniProt ID:	<u>Q8BXQ2</u>
Synonyms:	2510012P17Rik; 4930534E15Rik; CGI-06; NDAP; Ndap7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)