

## Product datasheet for **MP210867**

### Pitpnc1 Mouse qPCR Primer Pair (NM\_145823)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71795
Forward Sequence:	ATGTGACCGAGAAGGCGTGGAA
Reverse Sequence:	CGAGGTCTTTGGCTTCACTGTC
ACCN:	<u><a href="#">BC108351</a></u> , <u><a href="#">NM_145823</a></u> , <u><a href="#">NM_145823.1</a></u> , <u><a href="#">NM_145823.2</a></u> , <u><a href="#">BC028271</a></u> , <u><a href="#">BC082333</a></u>
UniProt ID:	<u><a href="#">Q8K4R4</a></u>
Synonyms:	1110020B03Rik; 5830436L09Rik; AI662802; AI851387; C330017I21Rik; Dnr411; rdgB-beta; RDGGB; RDGGB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)