

Product datasheet for **MP210875**

Pip4k2b Mouse qPCR Primer Pair (NM_054051)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 108083 |
| Forward Sequence: | CATGGTGTCCGAAACCTCAGG |
| Reverse Sequence: | GTCATAGGTGGTGAGGAAGCGA |
| ACCN: | BC047282 , NM_054051 , NM_054051.1 |
| UniProt ID: | Q80XI4 |
| Synonyms: | A1848124; c11; PI5P4Kbeta; Pip5k2b |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)