

Product datasheet for **MP210878**

Plcxd3 Mouse qPCR Primer Pair (NM_177355)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	239318
Forward Sequence:	CCCTCACCAATTTAGCCATCCC
Reverse Sequence:	CCGCATCAACTTCTTGGCTACAG
ACCN:	<u>NM_177355</u> , <u>NM_177355.1</u> , <u>NM_177355.2</u> , <u>NM_177355.3</u> , <u>BC137637</u> , <u>BC137639</u> , <u>BE945500</u>
UniProt ID:	<u>Q8BLJ3</u>
Synonyms:	B130016O10Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)