

Product datasheet for **MP210913**

Pacs2 Mouse qPCR Primer Pair (NM_001081170)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	217893
Forward Sequence:	CTCTGTGGACTACCGCTACAAC
Reverse Sequence:	GAGCACAGTTGGCTCCTGAGAT
ACCN:	<u>NM_001081170</u> , <u>NM_001081170.1</u> , <u>BC026674</u> , <u>BC043302</u> , <u>BC051499</u> , <u>BC059076</u>
UniProt ID:	<u>Q3V3Q7</u>
Synonyms:	6720425G15Rik; AW208651; AW489976; E230011J22; mKIAA0602; Pacs1l
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

