

## Product datasheet for **MP210915**

### **Pgm1 Mouse qPCR Primer Pair (NM\_028132)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	72157
<b>Forward Sequence:</b>	CAACGCACTGAAGGAGCTACTC
<b>Reverse Sequence:</b>	GGCACCAAGTTCTTCACAGAGG
<b>ACCN:</b>	<a href="#">BC008527</a> , <a href="#">BC086490</a> , <a href="#">NM_028132</a> , <a href="#">NM_028132.1</a> , <a href="#">NM_028132.2</a> , <a href="#">NM_028132.3</a> , <a href="#">BC005617</a> , <a href="#">BC036290</a> , <a href="#">BC055713</a> , <a href="#">BC066828</a> , <a href="#">BC080801</a>
<b>UniProt ID:</b>	<a href="#">Q9D0F9</a>
<b>Synonyms:</b>	2610020G18Rik; AA407108; AI098105; Pgm-2; Pgm1
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)