

Product datasheet for **MP210916**

Pgm2l1 Mouse qPCR Primer Pair (NM_027629)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70974
Forward Sequence:	GCTGTGGCAGAACTTCAGGAGA
Reverse Sequence:	AGAAGAGACAGTGGTGGCTAGC
ACCN:	BC032150 , BC076571 , NM_027629 , NM_027629.1 , NM_027629.2 , NM_027629.3 , BB633030 , BC020351 , BG071667
UniProt ID:	Q8CAA7
Synonyms:	4931406N15Rik; AI553438; BM32A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)