

Product datasheet for **MP210927**

Pik3ip1 Mouse qPCR Primer Pair (NM_178149)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216505
Forward Sequence:	GGTGCTACATCAGCAGCGAGAC
Reverse Sequence:	GGAACACCTGTGCCTCTTTGTC
ACCN:	<u>BC055920</u> , <u>BC079641</u> , <u>NM_178149</u> , <u>NM_178149.1</u> , <u>NM_178149.2</u> , <u>NM_178149.3</u> , <u>NM_178149.4</u> , <u>BC028479</u> , <u>BC053972</u>
UniProt ID:	<u>Q7TMJ8</u>
Synonyms:	1500004A08Rik; 5830455E04Rik; Crkd; Hgfl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

