

Product datasheet for **MP210943**

Pla2g4a Mouse qPCR Primer Pair (NM_008869)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18783
Forward Sequence:	GATGAGGCTCAAGGACCCAAAG
Reverse Sequence:	GAATAAAGCCGAGTCGCTCACC
ACCN:	<u>BC003816, NM_008869, NM_008869.1, NM_008869.2, NM_008869.3, NM_008869.4, BQ033204</u>
UniProt ID:	<u>P47713</u>
Synonyms:	cP; cPL; cPLA2; cPLA2-alpha; cPLA2alpha; Pla; Pla2g4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)