

Product datasheet for **MP210945**

Pla2g4c Mouse qPCR Primer Pair (NM_001004762)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	232889
Forward Sequence:	GGAGCAGTTTCTCGACCATCTG
Reverse Sequence:	CCACTGCCATTTGGAGACACAC
ACCN:	BC117808 , BC117809 , BC129974 , NM_001004762 , NM_001004762.1 , NM_001004762.2 , NM_001004762.3 , BC138548 , BC054740 , BC145767
UniProt ID:	Q64GA5
Synonyms:	D7Ert445e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)