

Product datasheet for **MP210947**

Pla2g4e Mouse qPCR Primer Pair (NM_177845)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	329502
Forward Sequence:	CAATGGCGAGAGGACCTCCAA
Reverse Sequence:	GCTGATCTGACAGCTTGCAATTCG
ACCN:	<u>NM_177845</u> , <u>NM_177845.1</u> , <u>NM_177845.2</u> , <u>NM_177845.3</u> , <u>NM_177845.4</u> , <u>BC137899</u> , <u>BC060600</u> , <u>BC066085</u> , <u>BC145020</u>
UniProt ID:	<u>Q50L42</u>
Synonyms:	2310026J01Rik; C230096D22; Pla2epsilon
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)