

## Product datasheet for **MP210948**

### Pla2g4f Mouse qPCR Primer Pair (NM\_001024145)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	271844
Forward Sequence:	TTCTGGATGCCGTGACCTACCT
Reverse Sequence:	GCAGACTCGCTCTGAGGCATAT
ACCN:	<a href="#">NM_001024145</a> , <a href="#">NM_001024145.1</a> , <a href="#">NM_001024145.2</a> , <a href="#">BC145426</a> , <a href="#">BC039947</a> , <a href="#">BC046400</a> , <a href="#">BC145427</a> , <a href="#">BC145428</a> , <a href="#">BC150771</a> , <a href="#">NM_001024145.3</a>
UniProt ID:	<a href="#">Q50L41</a>
Synonyms:	4732472I07Rik; Gm647; Pla2zeta
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)