

Product datasheet for **MP210950**

Pla2g6 Mouse qPCR Primer Pair (NM_016915)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53357
Forward Sequence:	AAGCCAGCGTTCATCCTGAGCT
Reverse Sequence:	TCGATGGCGATGAGAAGCTGGA
ACCN:	<u>BC003487</u> , <u>BC052845</u> , <u>BC057209</u> , <u>NM_016915</u> , <u>NM_016915.1</u> , <u>NM_016915.2</u> , <u>NM_016915.3</u> , <u>NM_016915.4</u> , <u>BC049778</u> , <u>BY099728</u>
UniProt ID:	<u>P97819</u>
Synonyms:	BB112799; iPLA(2)beta; iPLA2; iPLA2beta; PNPLA9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

