

## Product datasheet for **MP210957**

### Plcb2 Mouse qPCR Primer Pair (NM\_177568)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18796
Forward Sequence:	TGCTGGATGTCAGATGGTTGCC
Reverse Sequence:	CGTGCTTTAGGAGGTAGCCACT
ACCN:	<a href="#">NM_177568</a> , <a href="#">NM_177568.1</a> , <a href="#">NM_177568.2</a> , <a href="#">BB620972</a> , <a href="#">BB635904</a> , <a href="#">BC019788</a> , <a href="#">BC038048</a> , <a href="#">BC145249</a> , <a href="#">BC156752</a> , <a href="#">BE991828</a>
UniProt ID:	<a href="#">A3KGF7</a>
Synonyms:	AI550384; B230205M18Rik; B230399N12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)