

## Product datasheet for **MP210973**

### **Pld1 Mouse qPCR Primer Pair (NM\_008875)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	18805
<b>Forward Sequence:</b>	CTGGACTGCATCCTCAAACGGA
<b>Reverse Sequence:</b>	TCATCAACGTCCGCTTGCTGTAC
<b>ACCN:</b>	<u><a href="#">BC068144</a></u> , <u><a href="#">NM_008875</a></u> , <u><a href="#">NM_008875.1</a></u> , <u><a href="#">NM_008875.2</a></u> , <u><a href="#">NM_008875.3</a></u> , <u><a href="#">NM_008875.4</a></u> , <u><a href="#">BC026625</a></u>
<b>Synonyms:</b>	AA536939; C85393; mPLD1; Pld1a; Pld1b
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
<b>Storage:</b>	Store at -20 °C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)