

Product datasheet for **MP210983**

Ppcdc Mouse qPCR Primer Pair (NM_176831)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66812
Forward Sequence:	TCTGAAGCTGCCTCTCCTGGTA
Reverse Sequence:	TCATCAGCGTCACTGTAGAGGG
ACCN:	BC052928 , BC089351 , NM_176831 , NM_176831.1 , NM_176831.2 , NM_176831.3 , NM_176831.4 , BC004779 , BU053206 , BU694957 , NM_176831.5
UniProt ID:	Q8BZB2
Synonyms:	1810057I13Rik; 8430432M10Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)