

Product datasheet for **MP211003**

Phyhdl Mouse qPCR Primer Pair (NM_172267)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	227696
Forward Sequence:	TGTGCCTTTGCCAGTGAGGAGA
Reverse Sequence:	CTCCATGAGGTGGACAGTGTAG
ACCN:	<u>BC132271</u> , <u>BC132273</u> , <u>NM_172267</u> , <u>NM_172267.1</u> , <u>NM_172267.2</u> , <u>NM_172267.3</u> , <u>NM_172267.4</u> , <u>BC039982</u> , <u>BC145356</u> , <u>BC171925</u> , <u>BE945679</u> , <u>BQ829083</u> , <u>BU701903</u> , <u>BY196857</u> , <u>BY344206</u>
UniProt ID:	<u>Q9DB26</u>
Synonyms:	MGC49146
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

