

Product datasheet for **MP211015**

Pir Mouse qPCR Primer Pair (NM_027153)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69656
Forward Sequence:	GACGCTCAGAAAAGATGGTGGC
Reverse Sequence:	GGTGTGCGAGTGTAACCTTCG
ACCN:	BC024062 , NM_027153 , NM_027153.1 , NM_027153.2 , NM_027153.3 , BY504889 , BY718294 , BY720928
UniProt ID:	Q9D711
Synonyms:	2310042L19Rik; Pirnl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

