

Product datasheet for **MP211019**

Pttglip Mouse qPCR Primer Pair (NM_145925)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108705
Forward Sequence:	ACAAGGCGTGTATGGACTACCC
Reverse Sequence:	CAGGACCGACATGGTGATGATC
ACCN:	<u>BC023961</u> , <u>BC025533</u> , <u>BC029144</u> , <u>BC032184</u> , <u>NM_145925</u> , <u>NM_145925.1</u> , <u>NM_145925.2</u> , <u>NM_145925.3</u>
UniProt ID:	<u>Q8R143</u>
Synonyms:	1810010L20Rik; AI31431I; AU018448; C79540
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

