

Product datasheet for **MP211039**

Plat Mouse qPCR Primer Pair (NM_008872)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18791
Forward Sequence:	GTTACACAGCGTGGAGGACCAA
Reverse Sequence:	CAGTCAGCTTTTCGGTCCTCA
ACCN:	<u>BC011256</u> , <u>BC057967</u> , <u>BC061508</u> , <u>NM_008872</u> , <u>NM_008872.1</u> , <u>NM_008872.2</u> , <u>NM_008872.3</u>
UniProt ID:	<u>P11214</u>
Synonyms:	AU020998; AW212668; D8Ertd2; D8Ertd2e; t; t-; tPA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)